

DESCRIPTION

Transfer Bath Benches are designed to assist the user get into and out of the bathtub, it provides a sturdy base to sit on before sliding into position. The bench extends over the bath rim. Two legs with suction feet set inside of the bath prevents movement. The backrest and handrail attachments provide support for the user while transferring and bathing.

FEATURES

- Lightweight aluminium frame
- Height adjustable
- Plastic seat and back rest with non-slip textured surface
- Drainage holes in the seat reduces build up of water
- Reversible back rest to allow access from either side
- Handrail for support
- Wide base with suction cups tips for extra stability
- Minor assembly required
- Conditional warranty applies

SPECIFICATIONS

Seat height:	50.5 – 65.5cm
Maximum clearance:	55cm
Backrest height:	35cm
Armrest height:	29cm
Seating area:	68 x 41cm
Overall dimensions:	78.5 x 56 x 99.5cm (at highest setting)
Weight capacity:	200kg
Unit weight:	6kg

200kg Heavy Duty Transfer Bath Bench



SETTING UP

The legs with suction feet should be attached on the same side as the side rail, this side will be positioned closest to the wall inside the bath.

Determine which direction the transfer bench will face, the back rest can be attached on either side.

On each leg, press the spring button and slide the leg piece to the desired height.

Adjust each leg so the bench sits level. Adjust the chair to an appropriate seating height to allow the user to sit and stand with ease.

IMPORTANT: Ensure that the spring button has fully engaged the notch and the leg does not move.

HOW TO USE








Getting in the bathtub

1. To sit on the bench, move backwards until you can feel the seat behind your legs, then slowly lower yourself down and slide backwards slightly.
2. Once seated, turn (without twisting), raise and lower one leg at a time into the bath. Use the side rail for support as much as possible.
3. Slowly move along the bench as you move your legs in.

Getting out of the bathtub

1. When getting off the bench, ensure your feet are placed firmly on the floor and push up using your hands for support to stand leading with your head and shoulders.
2. Consider use of a non-slip mat in the bath under feet.

PRECAUTIONS AND SAFETY

-  Before use, check that the bench is in good working condition.
-  Be sure your weight is within the weight capacity of the unit and the product is suitable for your needs and abilities. If you are uncertain, seek the advice of a suitably qualified therapist or carer.
-  With all products, care must be taken to ensure stability and all feet are positioned firmly on the surface / floor.
-  Avoid tilting or leaning over too far.
-  DO NOT stand on the bench.
-  Use caution when having the bench at higher settings as it may become unstable, there is a greater risk of tipping backwards if the user sits down heavily.
-  Ensure that items required for bathing are easy to reach when seated.

MAINTENANCE AND CLEANING

- Check condition regularly for wear and tear.
- DO NOT use the equipment if it shows any sign of damage.
- Regularly clean the bench with a mild disinfectant.
- DO NOT use abrasive pads, cleaning agents or solvents.
- To prevent rust, dry product after use and before storing.
- Store in a dry place when not in use.