

40cm Portable Travel Transporter Wheelchair

FEATURES

- Powder coated aluminium frame
- Attendant locking hand brakes
- Seat belt
- Solid PVC castors
- Flip up armrest
- Black nylon upholstery with seat belt
- Folding aluminium footrest
- Conditional warranty applies

SPECIFICATIONS

PA112

Seat Width:	41cm
Seat Height:	45cm
Seat Depth:	38cm
Backrest Height:	36cm
Armrest Height:	28.5cm
Overall Width:	53cm
Folded Width:	37cm
Overall Length:	89cm
Overall Height:	83.5cm
Backrest folded:	58cm
Front Wheel Size:	12.5 x 2.5cm
Rear Wheel Size:	19.5 x 2.5cm
Weight Capacity:	100kg
Net Weight:	8kg



HANDLING AND SAFETY TIPS

Warning: Wheelchairs are heavy items – ensure proper lifting techniques to avoid injury. Ask for assistance whenever possible.

When lifting a wheelchair with the patient on it, do not hold on to the handles if the chair's backrest is the foldable kind as it may weaken / damage the folding hinge and bracket.

Be sure the user's weight is within the load capacity of the wheelchair and it is suitable for their needs and abilities. Seek the advice of a suitably trained specialist if uncertain.

Before using the wheelchair, always check to ensure product is in good working condition.

Regularly check brakes are in working order.

Ensure both brakes are applied before getting in or out of the wheelchair.

DO NOT stand on the footplates when getting in or out of the wheelchair.

DO NOT overload the wheelchair or use the wheelchair with the back folded or hang objects on the backrest – this can alter the centre of gravity and could cause the chair to tip backwards

Stability

The wheelchair is designed to be safe and stable if the centre of gravity is maintained while going about normal activities. Care must be taken while performing activities that shift the centre of gravity in the chair. Never reach so far that require you to shift in the seat or bend over out of the chair.

Ramps and Inclines

Users should be aware of their own capabilities, strength and limitations before attempting to go up a ramp. Correct techniques on ramp negotiation should be demonstrated by a qualified professional. Avoid changing direction when going down an incline as this could result in instability or tipping over. Seek assistance where possible

Caregivers should navigate downward slopes with the wheelchair turned to face up slope and back the chair downwards gradually.



1. Attendant hand grip
2. Attendant locking handbrake
3. Fold down backrest
4. Backrest fold down latch
5. Flip-up arm rests
6. Padded nylon seat*
7. 8" Rear castor with brake
8. 5" Front swivel castor
9. Folding footrest

**Photos for illustrative purposes only.*

Additional seat and back support cushion shown not included.

UNFOLDING the Wheelchair

Hold the attendant handgrips with both hands to separate the frame slightly. Then push down on the sides of the seat (6) until the frame is fully extended and seat locks into the seat guides.

FOLDING the Wheelchair

Fold the wheelchair by pulling up the centre of the seat (6), bring the sides together and swing footrests up; fold down the backrest by pressing the release latch (4).

BRAKE SYSTEM

Important: Ensure both brakes are applied before getting in or out of the wheelchair.

Attendant hand brake (2)

Squeeze brake handles (2) upward to apply brakes.

Push down brake handles (2) to lock brakes.

Squeeze handles (2) upward to release brake lock.

Note: When not in use, DO NOT leave hand brakes engaged as this may cause permanent indents on the tyre.

CARE AND MAINTENANCE

- Check condition regularly for wear and tear.
- Do not use if product shows signs of damage or corrosion.
- Always check bolts and nuts are securely fastened.
- Clean using a slightly damp cloth and wipe dry.
- Lubricate moving parts. We recommend a silicone based lubricant which will also assist in water / moisture displacement and can help reduce rust / corrosion.
- It is recommended that the wheelchair be serviced at least once a year by a mobility aid specialist.