

Wheelchair Contoured Seat Cushion

Product Code: BE1100 (suit a 16" wide seat)
BE1110 (suit a 18" wide seat)
BE1120 (suit a 20" wide seat)

- The Wheelchair Seat Cushion helps to stabilise and support the body.
- It is available in 3 seat sizes, 16" (40cm), 18" (45cm) and 20" (50cm).
- Its moulded contoured shape promotes comfort and proper body positioning.
- Made of a high quality PU foam.
- The cushion has a velcro hook and loop system which allows you to secure the cushion to a wheelchair seat.
- Is supplied with a removable fabric cover which is easy to clean, waterproof, anti-bacterial and fire retardant.



Specifications	BE1100	BE1110	BE1120
Width:	41cm	46cm	50cm
Depth:	40cm	40cm	40cm
Height:	6cm	6cm	6cm
Weight:	700gm	800gm	900gm



Pharmaquip Healthcare - 1/7 Monte Khoury Drive Queensland Australia 4129
PH:(07) 3806 4274 - www.pquip.com.au - sales@pquip.com.au