

Wheelchair Seat Cushion

Product Code: BE1160

- A general purpose seat cushion made from PU foam.
- It helps promote proper body positioning and provides a comfortable seating option.
- Supplied with a removable fabric cover which is easy to clean, waterproof , anti-bacterial and fire retardant.

Specifications

Width:	45cm
Depth:	40cm
Height:	6cm
Weight:	600gm



Pharmaquip Healthcare - 1/7 Monte Khourey Drive Queensland Australia 4129
PH:(07) 3806 4274 - www.pquip.com.au - sales@pquip.com.au