

Wheelchair Lumbar Cushion

Product Code: BE1400 (suit a 16" wide seat)
BE1410 (suit a 18" wide seat)
BE1420 (suit a 20" wide seat)

- Wheelchair Lumbar Cushions aid in improving your seating posture.
- The foam support conforms to the shape of the lumbar region to provide back support and comfort.
- The cushion has locking straps to secure it to the wheelchair backrest and a semi-rigid backing board .
- Is supplied with a fabric cover which is easy to clean, waterproof , anti-bacterial and fire retardant.



Specifications	BE1400	BE1410	BE1420
Width:	40cm	45cm	50cm
Height:	41cm	41cm	41cm
Thickness:	2cm-6cm-2cm	2cm-6cm-2cm	2cm-7cm-2cm
Weight:	800gm	900gm	1kg



Pharmaquip Healthcare - 1/7 Monte Khoury Drive Queensland Australia 4129
PH:(07) 3806 4274 - www.pquip.com.au - sales@pquip.com.au