

Description: **Crutch Underarm Pads Yellow**
Product Code: RACR102

- Designed to suit our underarm crutches, these pads are made from a high density natural sponge rubber providing complete comfort and avoiding underarm pressure.
- The interior is a corrugated cell construction for added comfort, anti-slip, and longer wear.
- Available in yellow only.
- Sold individually.

Specifications

Overall Length: 19cm



Description: **Crutch Hand Grips Yellow**
Product Code: RACR103

- Designed to suit our underarm crutches, these hand grips are soft and comfortable and will reduce wrist tension, finger cramp and hand blisters.
- Available in yellow only.
- Sold individually.

Specifications

Overall Length: 10cm

Overall Diameter: 4cm



Description: **22mm Crutch Tip**
Product Code: RACR100 - Grey
RACR101 - Yellow

- These replacement tips will suit our underarm crutches.
- They fit a shaft which is 22mm in diameter.
- Available in yellow and grey.
- Sold individually.

Specifications

Tip Diameter at Base: 4.3cm

Tip Height: 5.6cm



Crutch Spare Parts

Description: **22mm Heavy Duty Crutch Tip Grey**
Product Code: **RTIP101**

- These are a heavy duty replacement tips that our suit our underarm and forearm crutches.
- They fit a shaft which is 22mm in diameter.
- Available in grey only.
- Sold individually.

Specifications

Tip Diameter at Base: 5cm
Tip Height: 7.3cm



Description: **22mm Forearm Crutch Tip Grey**
Product Code: **RTIP103**

- These replacement tips will suit our forearm crutches.
- They fit a shaft which is 22mm in diameter.
- Available in grey only.
- Sold individually.

Specifications

Tip Diameter at Base: 4cm
Tip Height: 5.9cm



Description: **8mm Brass Button and Spring**
Product Code: **SCR300**

This 8mm diameter brass button and spring will suit the following series crutches:
RCR101 to RCR103, RCR301 to RCR303



Description: **6.5mm Brass Button and Double Spring**
Product Code: **SCR310**

This part is supplied with two 6.5mm brass buttons and one double spring and will suit the following series crutches:
RCR501 to RCR503



Crutch Spare Parts

Description: **Crutch Plastic Hand Grip**
Product Code: **SCR100**

Includes the plastic hand grip, nut, bolt and washers, to suit our range of Underarm Crutches:
RCR101 to RCR103 Series

Handgrip Length: 10.5cm
Overall Diameter: 2.5cm



Description: **Forearm Cuff and Stem - Vinyl**
Product Code: **SCR200**

Replacement cuff and stem for Forearm Crutches:
RCR301 to RCR303 Series

Note: The forearm cuffs for the RCR301-303 series and RCR501-503 series are not interchangeable.



Description: **Forearm Cuff and Stem - Plastic**
Product Code: **SCR210**

Replacement cuff and stem for Forearm Crutches Ergonomic
Handgrip:
RCR501 to RCR503 Series

Note: The forearm cuffs for the RCR301-303 series and RCR501-503 series are not interchangeable.

