

Toilet Surround

Product Code: PQ103

The Toilet Surround is a sturdy powder coated, height adjustable, steel frame, designed to assist patients in lowering themselves onto the toilet, and raising from it. The frame sits around the toilet and can easily be removed when not in use.

Features

- Powder coated white finish
- Height adjustable
- Doesn't require fixing to toilet
- 12 months warranty (conditions apply)

Specifications

Handle Height:	63cm to 76cm
Width Between Arms:	48cm
Width at Base:	53cm
Depth at Base:	44cm
Weight Capacity:	114kg
Net Weight:	3.4kg



Warranty

- Pharmaquip Healthcare offers a one year conditional warranty for domestic use of this product.
- If this product is used for commercial purposes (for example, hire), the warranty period is limited to 6 months.
- Pharmaquip Healthcare warrants its products to be free from defects in material and workmanship.
- Pharmaquip Healthcare will repair or replace any part, provided the product is delivered, pre-paid, to a location authorised by Pharmaquip Healthcare within the prescribed period.
- This warranty does not include parts subject to wear and tear.
- This warranty does not cover damage caused by misuse, neglect, accident or improper use.
- Warranties issued by Pharmaquip Healthcare will be limited to repair or replacement of the defective product.
- Pharmaquip Healthcare will not be liable for any consequential or incidental damages whatsoever.
- This product is intended to be used as a toilet aid. Pharmaquip Healthcare will not be liable for any other improper use.