

# Bath Seat

Product Code: RBN201 without Back  
RBN202 with Back



Bath Seats are compact and durable and are designed to be placed in the bath or shower to provide safety and comfort while bathing. Made of aluminium, with a textured plastic blow moulded seat, they are height adjustable and have slip-resistant rubber tips.

## Features

- Height adjustable in 25.4mm (1") increments
- Suction cup tips on all four legs of the bath bench provide extra stability and safety
- High density curved polyethylene seat with handgrips
- Durable aluminium frame is corrosion resistant
- Drainage holes in seat help reduce the build up of water
- Angled legs distribute weight over large area
- Assembly required
- 12 months warranty (conditions apply)

## Specifications

Seat Height:	35cm to 52.5cm
Seat Width:	50cm
Seat Depth:	31cm
Backrest Height:	68cm to 85cm (RBN202)
Width at Base:	42cm
Depth at Base:	44cm
Weight Capacity:	114kg
Net Weight:	2kg (RBN201) 3.4kg (RBN202)

## Spare Parts



RTIP104  
Replacement tip

