

## Bath Stool

Product Code: RBN203

Ideal for narrow bath tubs or shower areas. This height adjustable stool is a lightweight, aluminium frame with a plastic seat and anti-slip rubber tips.

### Features

- Height adjustable in 25.4mm (1") increments
- Suction cup tips on all four legs provide extra stability and safety
- Corrosion resistant aluminium frame
- Supplied fully assembled.
- 12 months warranty (conditions apply)

### Specifications

Seat Height: 39cm to 49cm  
Seat Diameter: 31cm  
Width at Base: 32.5cm  
Depth at Base: 32.5cm  
Weight Capacity: 114kg  
Net Weight: 2kg



## Care and Maintenance

- Check rubber tips for tears or cracks, replace as necessary.
- Check the leg extensions are securely locked in.
- Ensure all screws, nuts, bolts etc are securely tightened.
- Regularly clean with warm soapy water. For tougher stains use a light cream cleanser.

## Spare Parts

Replacement tip: RTIP104

## Warranty

- Pharmaquip Healthcare offers a one year conditional warranty for domestic use of this product.
- If this product is used for commercial purposes (for example, hire), the warranty period is limited to 6 months.
- Pharmaquip Healthcare warrants its products to be free from defects in material and workmanship.
- Pharmaquip Healthcare will repair or replace any part, provided the product is delivered, pre-paid, to a location authorised by Pharmaquip Healthcare within the prescribed period.
- This warranty does not include parts subject to wear and tear.
- This warranty does not cover damage caused by misuse, neglect, accident or improper use.
- Warranties issued by Pharmaquip Healthcare will be limited to repair or replacement of the defective product.
- Pharmaquip Healthcare will not be liable for any consequential or incidental damages whatsoever.
- This product is intended to be an aid in bathing in the bath or shower. Pharmaquip Healthcare will not be liable for any other improper use.