

Grab Bars

Product Code: RGB201 - 12" (300mm)
 RGB202 - 16" (406mm)
 RGB203 - 18" (457mm)
 RGB204 - 24" (610mm)
 RGB205 - 32" (813mm)

The RGB series of grab bars are a textured grab bar which helps provide safety when getting into or out of the bathtub, or assistance when using the toilet.

Features

- Made of a heavy gauge chrome plated 1" steel tubing
- 76mm Diameter flanges
- Grab Bars sit 38mm out from the wall
- Knurled chrome plated finish reduces the danger of hands slipping
- Mounts vertically, horizontally or diagonally
- Easy to install, however professional installation is highly recommended
- Mounting hardware included
- 12 months warranty (conditions apply)



Grab measurements are based on the distance between the centre of the flanges



Warranty

- Pharmaquip Healthcare offers a one year conditional warranty for domestic use of this product.
- If this product is used for commercial purposes (for example, hire), the warranty period is limited to 6 months.
- Pharmaquip Healthcare warrants its products to be free from defects in material and workmanship.
- Pharmaquip Healthcare will repair or replace any part, provided the product is delivered, pre-paid, to a location authorised by Pharmaquip Healthcare within the prescribed period.
- This warranty does not include parts subject to wear and tear.
- This warranty does not cover damage caused by misuse, neglect, accident or improper use.
- Warranties issued by Pharmaquip Healthcare will be limited to repair or replacement of the defective product.
- Pharmaquip Healthcare will not be liable for any consequential or incidental damages whatsoever.
- This product is intended to be an aid in bathing in the bath or shower. Pharmaquip Healthcare will not be liable for any other improper use.