

# Wheelchair Amputee Support

Product Code: SWH465L-Left  
SWH465R-Right

- The Amputee Support provides a comfortable cushioning support with a soft moulded foam pad.
- When placed on the wheelchair, the distance between the edge of the seat and the support can be adjusted from 8cm to 10.50cm to 13cm.
- The height of the support can also be adjusted to suit individual needs.
- This support is designed to suit our Lightweight Wheelchair, model PA210.



Pharmaquip Healthcare - 1/7 Monte Khoury Drive Loganholme Queensland Australia 4129  
PH:(07) 3806 4274 - [www.pquip.com.au](http://www.pquip.com.au) - [sales@pquip.com.au](mailto:sales@pquip.com.au)