

Handy Grip Reacher

Product Code: YDL310 - 60cm Small
YDL309 - 75cm Medium
YDL311 - 90cm Large

This lightweight reacher has an easy to use trigger handle making it ideal for use around the home.

Available in 3 sizes.

Features

- Magnetic head which allows user to pick up small metal objects such as pins
- Hook on the end of the reacher which aids in dressing
- Palm lever which helps clamp the jaws for those who have poor hand strength
- Handy clip to attach to a rollator or wheelchair
- Rotating head to reach awkward items



Magnetic head for small metal objects



Hook can be used as a dressing aid



Palm lever can be used to assist in clamping the gripping jaws



Clip to attach the reacher to a wheelchair or rollator



The gripping jaws can be rotated to aid in reaching into awkward places