

# Cane Tips

## Cane Tips

Product Code: RTIP107BL Black  
RTIP107BR Bronze  
RTIP107GR Grey

These replacement cane tips will fit a standard cane that has a 19mm diameter shaft (on the lower shaft). Sold individually.



## Large Base Quad Cane Tips

Product Code: RTIP109BL Black  
RTIP109GR Grey

These replacement cane tips will fit a large base quad cane that has a 16mm diameter shaft. Sold individually.



## Small Base Quad Cane Tips

Product Code: RTIP110BL Black  
RTIP110GR Grey

These replacement cane tips will fit a small base quad cane that has a 13mm diameter shaft. Sold individually.



## Cane Seat Tips

Product Code: RTIP113

These replacement cane tips will fit our tripod cane seat that has a 19mm diameter shaft. Sold individually.



## Specifications

Part No.	Colour	Suit Shaft Diameter	Tip Diameter At Base	Tip Height
RTIP107BL	Black	19mm	3.80cm	4.70cm
RTIP107BR	Bronze	19mm	3.80cm	4.70cm
RTIP107GR	Grey	19mm	3.80cm	4.70cm
RTIP109BL	Black	16mm	3.10cm	3.60cm
RTIP109GR	Grey	16mm	3.10cm	3.60cm
RTIP110BL	Black	13mm	2.80cm	2.90cm
RTIP110GR	Grey	13mm	2.80cm	2.90cm
RTIP113	Black	19mm	2.85cm	3.30cm

