

19mm Cane Tips

Product Code: RTIP107BL Black
 RTIP107BR Bronze
 RTIP107GR Grey



These replacement cane tips will fit a standard cane that has a 19mm diameter shaft (on the lower shaft).
 Sold individually.

16mm Quad Cane Tips

Product Code: RTIP109BL Black
 RTIP109GR Grey



These replacement cane tips will fit a large base quad cane that has a 16mm diameter shaft.
 Sold individually.

13mm Quad Cane Tips

Product Code: RTIP110BL Black
 RTIP110GR Grey



These replacement cane tips will fit a small base quad cane that has a 13mm diameter shaft.
 Sold individually.

Specifications

Part No.	Colour	Suit Shaft Diameter	Tip Diameter At Base	Tip Height
RTIP107BL	Black	19mm	3.80cm	4.70cm
RTIP107BR	Bronze	19mm	3.80cm	4.70cm
RTIP107GR	Grey	19mm	3.80cm	4.70cm
RTIP109BL	Black	16mm	3.10cm	3.60cm
RTIP109GR	Grey	16mm	3.10cm	3.60cm
RTIP110BL	Black	13mm	2.80cm	2.90cm
RTIP110GR	Grey	13mm	2.80cm	2.90cm