

Rollator Spare Parts

Bags, Baskets and Trays

Description: **Vinyl Bag (Under Seat)**
Part No: SRO400
Width: 31cm
Depth: 24cm
Height: 14cm



Suit Models: PA323A:PA323B; PA323C; PA323D

Description: **Vinyl Bag (3 Wheel)**
Part No: SRO410
Width: 32cm
Depth: 20cm
Height: 34cm



Suit Models: PA323H

Description: **Wire Basket (Under Seat)**
Part No: SRO420
Width: 37cm
Depth: 23cm
Height: 14.50cm



Suit Models: PA323A:PA323B; PA323C; PA323D

Description: **Wire Basket (Front)**
Part No: SRO430
Width: 40cm
Depth: 30cm
Height: 28cm



Suit Models: PA323E

Description: **Wire Basket (Large)**
Part No: SRO440
Width: 52cm
Depth: 21cm
Height: 18cm



Suit Models: RRO801

Description: **Rollator Seat Tray**
Part No: SRO460
Width: 39.50cm
Depth: 26cm

Suit Models: PA323E



Rollator Spare Parts

Brake Assembly

Product Code	Description	Brake cable To suit	Suit Models
SRO101L	Complete Brake Assembly Left-PA323A,B,D	SRO123	PA323A, PA323B, PA323D
SRO101R	Complete Brake Assembly Right-PA323A,B,D	SRO123	PA323A, PA323B, PA323D
SRO102L	Complete Brake Assembly Left-PA323E	SRO122	PA323E
SRO102R	Complete Brake Assembly Right-PA323E	SRO122	PA323E
SRO103L	Complete Brake Assembly Left-PA323H	SRO123	PA323H
SRO103R	Complete Brake Assembly Right-PA323H	SRO123	PA323H
SRO108L	Complete Brake Assembly Left-RRO801	SRO121	RRO801
SRO108R	Complete Brake Assembly Right-RRO801	SRO121	RRO801



SRO101L/R



SRO102L/R



SRO103L/R



SRO108L/R



Rollator Spare Parts

Brake Cables


Description: **Brake Cable-88cm**
 Part No: SRO123
 Overall Length: 88cm
 Black Sheath: 80cm
 Suit Models: PA323A; PA323B; PA323D;



Description: **Brake Cable-110cm**
 Part No: SRO121
 Overall Length: 110cm
 Black Sheath: 98cm
 Suit Models: RRO801



Description: **Brake Cable-117cm**
 Part No: SRO122
 Overall Length: 113cm
 Black Sheath: 102cm-with centre tensioner
 Suit Models: PA323E



Description: **Brake Cable Ties**
 Part No: SRO128



Rollator Brake Levers

Description: **Brake Lever System
Left and Right**
 Part No: SRO130L-Left
 SRO130R-Right
 Suit Models: PA323A; PA323B;
 PA323D; PA323H



Description: **Rear Brake Lever**
 Part No: SRO136L-Left
 SRO136R-Right
 Suit Models: RRO801



Description: **Rear Brake**
 Part No: SRO133
 Suit Models: PA323E



Description: **Spring Loader**
 Part No: SRO134
 Suit Models: PA323C



Rollator Spare Parts

Castors



Round Profile



Flat Profile



Description:	5" Castor	6" Castor	8" x 1 1/4" Castor	8" x 2" Castor
Product Code:	SRO200	SRO201	SRO202	SRO203
Profile:	Round	Round	Flat	Round
Diameter:	123mm (5")	148mm (6")	194mm (8")	188mm (8")
Width:	25.4mm (1")	25.4mm (1")	34mm (1.1.4")	50mm (2")
Bearing:	608zz	608zz	608zz	608zz
Hub:	36mm	36mm	38mm	60mm
Suit Models:	PA323C (Front castor)	PA323B	PA323D	PA323E:PA323H

Castor Forks



Description: **Front Castor Fork - 5"/6" Wheel**
 Product Code: SRO220
 Includes: Castor Fork; Fork Screw and Bearing
 Suits Models: PA323A;PA323B;PA323C

Description: **Front Castor Fork - 8" Wheel**
 Product Code: SRO221
 Includes: Castor Fork; Fork Screw and Bearing
 Suits Models: PA323D



Description: **Front Castor Fork - 8" Wheel (801)**
 Product Code: SRO222
 Includes: Castor Fork; Fork Screw and Bearing
 Suits Models: RRO801

Description: **Front Castor Fork - 14 x 6 x 12**
 Product Code: SWH720
 Suits Models: PA323E; PA323H
 Chrome finish only



Rollator Spare Parts

Castor Bearings

Description: **Castor Bearing 608ZZ (22 x 7 x 8)**

Product Code: SRO210

Centre: 8mm
Diameter: 22mm
Width: 7mm



Suits Models: PA323A;PA323B; PA323C (front wheel); PA323D; PA323E; PA323H; RRO801

Castor Fork Bearings

Description: **Castor Fork Bearing 6001Z (28 x 8 x 12)**

Product Code: SWH601

Centre: 12mm
Diameter: 28mm
Width: 8mm



Suits Models: PA323A;PA323B; PA323C; PA323D; PA323E

Description: **Castor Fork Bearing 30/27 x 9 x 12**

Product Code: SWH606

Centre: 12mm
Diameter: 30mm/27mm
Width: 9mm



Suits Models: PA323H

Axle Dust Cap

Description: **Castor Axle Dust Cap**

Product Code: SRO238



Rollator Castor Fork Axle Bolts

Description:	Front Castor Axle Bolt- 65mm	Rear Castor Axle Bolt- 75mm	Front Castor Axle Bolt- 75mm	Rear Castor Axle Bolt- 114mm	Front Castor Axle Bolt- 78mm	Rear Castor Axle Bolt- 105mm
Product Code:	SRO230	SRO231	SRO232	SRO233	SRO234	SRO235
Length:	65mm	75mm	75mm	114mm	78mm	105mm
Diameter:	8mm	8mm	8mm	6mm	8mm	8mm
Type:	Hex Head	Hex Head	Dome Head	Hex head	Dome head	Dome Head
Suit Models:	PA323A; PA323B; PA323C; PA323D	PA323A; PA323B; PA323D	PA323E; PA323H	PA323C (rear wheel only)	RRO801	RRO801



Rollator Spare Parts

Handgrips

Description: **Orthopaedic Hand Grip**
Part No: SRO300L-Left
SRO300R-Right
Suit Models: Will suit any rollator handle 22mm
in diameter



Description: **Foam Handgrip**
Part No: SRO320
Suit Models: Will suit any rollator handle 22mm
in diameter



Height Adjustment Knobs

Description: **Height Adjustment Knob and
Screw**
Part No: SRO500
Suit Models: PA323A; PA323B; PA323C;
PA323D



Description: **Height Adjustment Lever and
Screw**
Part No: SRO510
Suit Models: PA323E; PA323H

Note: Two screws are supplied with the lever.
The 30mm screw suits the PA323H model and the
37mm screw suits the PA323E model



Rollator Spare Parts

Backrests

Description: **Pushbutton Detachable Backrest**
Part No: SRO700
Width: 49.5cm
Suit Models: PA323B; PA323C; PA323D
Colours: Red, Blue, Green and Silver



Description: **Pushbutton Detachable Backrest (small)**
Part No: SRO710
Width: 44.5cm
Suit Models: PA323A
Colours: Red, Blue, Green and Silver



Description: **Pushbutton Detachable Backrest (large)**
Part No: SRO720
Width: 62cm
Suit Models: RRO801
Colours: Blue only



Description: **Detachable Backrest Folding Bracket**
Part No: SRO740
Suit: SRO700 and SRO710 Backrests



Rollator Spare Parts

Seats

Description: **Padded Seat**
Part No: SRO600
Seat Width: 31cm
Seat Depth: 31cm
Suit Models: PA323A; PA323B; PA323C;
PA323D



Seat hinges not included

Description: **PU Seat**
Part No: SRO610
Seat Width: 38cm
Seat Depth: 16cm
Suit Models: PA323E
(PU = polyurethane)



Description: **Large Padded Seat**
Part No: SRO620
Seat Width: 46cm
Seat Depth: 33cm
Suit Models: RRO801



Seat hinges not included

Description: **Seat Hinge and Screws**
Part No: SRO630
Suit Models: SRO600 and SRO620 seats

