

Walking Frame Spare Parts

28.4mm Walker Tips

Product Code:
RTIP105BL- Black
RTIP105GR- Grey

These replacement tips will fit a walker frame that is 28.4mm in diameter on the extension leg.

Sold individually.

Suit Models: RBR601, RBR602, RWK101, RWK301, RWK303, RWK502, RWK503

Tip Diameter at Base: 4cm

Tip Height: 4cm



25.4mm Walker Tips

Product Code:
RTIP106BL- Black
RTIP106GR- Grey

These replacement tips will fit a walker frame that is 25.4mm in diameter on the extension leg.

Sold individually.

Suit Model: RWK402

Tip Diameter at Base: 4cm

Tip Height: 4cm



Forearm Pad

Product Code: SWF150

These forearm pads are designed to suit our PA65 Gutter Arm walker. The pad is vinyl coated with soft padding and a velcro pad to attach to the walking frame.

Sold individually.

Length: 23cm

Width: 16cm

Thickness: 2.5cm



Armrest Height Adjustment Knob

Product Code: SWF130

Adjustment knob for the armrest on our PA65 Gutter Arm Walker



Plastic End Caps

Product Code:
SWF110- 25.4mm
SWF120- 28.4mm

Plastic Ends Caps are designed to fit on the end of the extension leg tube. They are available in two sizes:

- SWF110 is a 25.4mm diameter tube and is suitable for our model RWK402.
- SWF120 is a 28.4mm diameter tube and suits all other models.

