

Wheelchair Elevating Legrests

Elevating Legrests provide the option to elevate the leg from an angled position to a horizontal position.

These Legrests have an adjustable calf pad to support the leg and the footplate can be extended to protect the foot.

Once the Legrest is attached to the chair, simply raise it to the desired level.

A release lever on the side of the Legrest will lower the leg down.

Available for either the left or right side.

Legrests from Pharmaquip Healthcare are designed to suit our range of wheelchairs and may not be compatible with other brands.



Description: **Wheelchair Elevating Legrest (A)**

Part No: SWH454L-Left

SWH454R-Right

Suit Models: PA160, PA190, PA195



Description: **Wheelchair Elevating Legrest (B)**

Part No: SWH456L-Left

SWH456R-Right

Suit Models: PA130



Description: **Wheelchair Elevating Legrest (C)**

Part No: SWH458L-Left

SWH458R-Right

Suit Models: PA200



Description: **Wheelchair Elevating Legrest (D)**

Part No: SWH459L-Left

SWH459R-Right

Suit Models: PA210



Pharmaquip Healthcare - 1/7 Monte Khoury Drive Loganholme Queensland Australia 4129

PH:(07) 3806 4274 - www.pquip.com.au - sales@pquip.com.au